

Viscotears® HA Plus

Preservative free eye drops containing 0.2% sodium hyaluronate and 0.2% hypromellose

10ml Multi-dose bottle

CAREFULLY READ THIS LEAFLET BEFORE FIRST USE

Viscotears® HA Plus is a sterile, ophthalmic solution with a physiological pH which provides lasting moisturisation and lubrication. It is based on sodium hyaluronate and hypromellose.

The combination of sodium hyaluronate and hypromellose wets and lubricates the ocular surface, forms a protective viscoelastic shield and re-establishes the tear film. Viscotears® HA Plus provides lasting relief from the symptoms of moderate to severe of dry eyes such as redness, burning, irritation and foreign-body sensation. These symptoms are often associated with contact lens wear, the use of computer or device screens, certain medical conditions and treatments, environmental factors such as air-conditioning or central heating and atmospheric pollution.

Thanks to the OSD closure system, which maintains the microbiological, chemical and physical characteristics, Viscotears® HA Plus does not contain any preservatives. The absence of preservatives allows this product to be used by contact lens wearers without the need to remove their lenses first.

INSTRUCTIONS FOR USE

- Tear off the tamper-evident seal before the first use.
- Open the bottle by removing the protective cap.
- Apply 1-2 drops in each eye by gently squeezing the bottle.
- After use, close the bottle: the contents should be used within 6 months from the date of first opening.
- Suitable for frequent use: Unless instructed otherwise by your doctor or pharmacist, use these drops 2-3 times per day or as often as is required. As dry eye is a condition which varies from person to person, the drops should be used as often as you think necessary, to relieve your dry eye.

As the bottle empties during the course of normal use, you may need to apply slightly more pressure when you squeeze the bottle to deliver the drops. This is due to the OSD closure system that is used to preserve the sterility of the preservative free solution.

INGREDIENTS

0.2% sodium hyaluronate, 0.2% hypromellose, boric acid, sodium tetraborate decahydrate, sodium chloride, W.F.I.

Free from preservatives and phosphates.

BAUSCH + LOMB

WARNINGS AND CAUTIONS

- The product is for single patient use.
- The product is only intended for ophthalmic use.
- Do not use this product if you are hypersensitive to any of the ingredients.
- Do not use the product if, at the time of first opening, the product appears damaged or the tamper evident seal is broken.
- During application, do not touch your eyes or any other surface with the tip of the bottle.
- If irritation occurs, stop using the product and speak to your doctor.
- Do not use after the expiry date stated on the carton and multi-dose bottle.
- Dispose of the container after use following local regulations.
- Keep out of reach of children.

COMMON SIDE EFFECTS

- The product can cause short term blurring of vision that clears after blinking. Do not drive or use machinery until your vision is clear.

If symptoms persist or worsen, contact your doctor.

PRESENTATION

10ml multi-dose bottle, preservative free.

STORAGE CONDITIONS

Do not store above 30°C.

Discard 6 months after first opening.



C.O.C. Farmaceutici S.r.l.
Via Modena, 15
40019 - Sant'Agata
Bolognese (Bo)
Italy

MEDICAL DEVICE

CE 0546

DISTRIBUTED BY:

Bausch & Lomb UK Ltd
106 London Road
Kingston upon Thames
KT2 6TN

SYMBOLS



Sterilised using aseptic processing techniques



Upper temperature limit



Do not use if packaging is damaged



Consult instructions for use



Caution: Consult the "instructions for use" for important cautionary information such as warnings and precautions



Use by date



Batch code



Manufacturer

®/™ are trademarks of Bausch & Lomb Incorporated or its affiliates

© Bausch & Lomb Incorporated

Date of last revision 03.01 of 07/2017 BLv1

BAUSCH + LOMB